

TOOTHPRINTS

1 - IDENTIFICATION

Distributor: Elevate Oral Care
Address: 346 Pike Road, Suite 5
City, State, Zip: West Palm Beach, FL 33411
Telephone: 877-866-9113
8:30am-5:00pm Eastern Time

Date Prepared: February 6, 2019

2 - COMPOSITION INFORMATION**Hazardous Ingredients**

This product contains no hazardous components as defined in the OSHA Hazard Communication Standard (29 CFR 1910.1200).

3 - PHYSICAL AND CHEMICAL PROPERTIES

Boiling Point: N/A
Vapor Pressure (mm Hg): N/A
Vapor Density (AIR = 1): N/A
Solubility in Water: Insoluble
Appearance and Odor: White, slight odor
Specific Gravity (H₂O = 1): Above 1.0
Melting Point: 140 °F – 144 °F (60 °C – 62 °C)

4 - FIRE AND EXPLOSION HAZARD DATA

Flash Point (Method Used): 380 ° ±5 °C
Extinguishing Media: Water spray, carbon dioxide, dry chemical, universal type foam
Special Fire Fighting Procedures: Fire Fighters should wear full protective clothing including a self-contained breathing apparatus.

5 - REACTIVITY DATA

Stability: Stable
Conditions to Avoid: None
Incompatibility (Materials to Avoid): None
Hazardous Decomposition or By-products: Will not occur
Hazardous Polymerization: Will not occur
Conditions to Avoid: None

6 - HEALTH HAZARD DATA

Routes of Entry:
Skin: None currently known
Eyes: None currently known
Inhalation: No emergency care anticipated
Ingestion: None known
Carcinogenicity - NTP: No
IARC Monographs: No OSHA Regulated Carcinogen: No

7 - EMERGENCY FIRST AID PROCEDURES

Eyes: Flush promptly with flowing water. If irritation develops, consult a physician.
Skin: Wash promptly with soap and water.
Inhalation: No emergency care anticipated.
Ingestion: No harmful effects expected.

8 - PRECAUTIONS FOR SAFE HANDLING & USE

Steps to be taken in case material is released or spilled: Collect for disposal
Waste Disposal Method: Disposal must be carried out in accordance with federal, state and local regulations.
Precautions to be taken in handling and storing: N/A
Other Precautions: N/A

9 - CONTROL MEASURES

Respiratory Protection: Not required
VENTILATION:
Local Exhaust: General room ventilation is satisfactory
Mechanical (General): N/A
Protective Gloves: Not required
Eye Protection: Safety glasses
Other Protective Clothing or Equipment: Eye bath
Work/Hygiene Practices: Avoid breathing dusts

10 - TRANSPORTATION INFORMATION

Not DOT regulated.

11 - SPECIAL INFORMATION

HMIS (Hazardous Material Identification System) Rating:
H1 F0 R0
[HMIS Hazard Index: 4 – Severe Hazard; 3 – Serious Hazard;
2 – Moderate Hazard; 1 – Slight Hazard; 0 – Minimum Hazard]